# YUWELL

## Mesh Nebulizer

### (Model:M1 02)

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v. Symbols regarding safety requirement of the machine and meaning concerned

Symbol	Meaning	Symbol	Meaning
	Type B applied part		Attention! Please consult document attached!
	Prohibit		Attention
	Fragile article		Lift up
	Required to be protected from rain		

vi. Common Fault Analysis and Troubleshooting

Fault	Cause Analysis	Troubleshooting
	Medication cup module is not assembled correctly	Restart after medication cup module is reassembled
Too small volume of spray	Liquid is used up or does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	Spray sheet is blocked	Clean medication cup according to operation instructions; change a new cup if the fault is not removed through cleaning

Fault	Cause Analysis	Troubleshooting
	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode
After power on, the power indicator LED flashes to off	Medication cup module is not assembled correctly	Restart after medication cup module is reassembled
	Medication cup has no liquid or liquid is used up	Inject liquid
	Liquid does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode
Power indicator LED is not on or nebulizer does not work	Reverse connection of battery anode and cathode	Reinstall battery according to operation instructions
	Too low battery power	Restart after changing a new battery
	Indicator LED (orange LED) is always on and battery power is insufficient	Restart after changing new batteries.
Power indicator LED is on but nebulizer does not work	The electrode of main unit or medication cup is dirty	Restart after cleaning electrode
	Spray sheet is dirty or seriously blocked	Clean medication cup according to operation instructions; change a new cup if the fault is not removed through cleaning

Fault	Cause Analysis	Troubleshooting
Nebulizer powers off automatically while it is operating	Medication cup module is not assembled well and therefore becomes loose in usage process	Restart after reassembling medication cup
	Medication cup has no liquid or liquid is used up	Inject liquid again
	Liquid does not contact spray sheet for longer than 10s	Incline the positive surface of nebulizer to the user and make liquid contact spray sheet
	Nebulizer is shaken or vibrated seriously when using it	Hold nebulizer by hands stably when using it
	Fault of medication cup	Change a new medication cup
Nebulizer cannot power off automatically	The electrode of main unit or medication cup is dirty	Power off nebulizer to clean electrode
	Fault of medication cup module	Change a new medication cup module
Liquid leakage	Medication cup module is damaged or sealing silica gel becomes aged	Change a new medication cup module

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JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO. •LTO.

Sales Center Address: Huanyuan East Road No.1. Xuzhuang Software Park, Nanjing, Jiangsu Province,

P.R.China.210000 Manufacturing Address: Yun Yang Industria1 Park. Dan Yang, Jiangsu Province, P.R. China, 212300

<http://www.yuwell.com>

Διανομή: [www.alfacare.gr](http://www.alfacare.gr)