EN | Caution: Please read before using

Consult a physician prior to use.

If contents comes in contact with eyes or skin, flush with water and contact a physician.

Never apply a Hot/Cold Pack directly to skin. This may cause burns.

Do not use on sensitive skin areas or on an open wound.

People with known or suspected circulatory problems should not use a Hot/Cold Pack except under direct supervision of a doctor.

Do not use a Hot/Cold Pack along with sports creams, balms, liniments or lotions.

Unattended use by children or incapacitated persons may be dangerous.

This is not a toy.

Do not use a Hot Pack on any part of the body where there is swelling, redness or immediately after injury. Cold therapy is recommended for these conditions.

For longer product life:

Use the Hot/Cold Pack for cold or for hot use, but not for both.

When placing the pack in the freezer, lay it flat, not folded.

Place the Hot/Cold Pack on top of you. Not the other way around.

COLD THERAPY:

Place gel pack into freezer for at least 1 to 2 hours.

Wrap Hot/Cold Pack in a cloth or towel, or use the special MoVeS Hot/Cold Pack Protection Sleeve.

Remove Hot/Cold Pack if it becomes uncomfortable.

Do not apply for more than 20 minutes at a time.

HOT THERAPY:

Microwave:

Place the Hot/Cold Pack flat in microwave oven.

Heat for 1 minute (750W).

Repeat with 30 seconds intervals until the desired temperature is achieved.

Hot Water:

Bring water to boil and remove from heat.

Place Hot/Cold Pack in hot water for 2 minutes or more until desired temperature is reached.

Remove Hot/Cold Pack from water with a pair of tongs.

Knead to distribute heat evenly.

If the Hot/Cold Pack is too hot, let cool before using.

Wrap Hot/Cold Pack in a cloth or towel, or use the special MoVeS Hot/Cold Pack Protection Sleeve.

Remove Hot/Cold Pack if it becomes uncomfortable.

Do not apply for more than 20 minutes at a time.

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