

AeroTrach[®] Plus^{*}

Anti-Static Valved Holding Chamber



Consistent and comfortable medication delivery for your tracheostomy patients

Effective Drug Delivery Right Out of the Package


The **AeroTrach Plus^{*}** Anti-Static Valved Holding Chamber (VHC) can be used directly out-of-package without pre-treatment – saving valuable time in hospital environments. The **AeroTrach Plus^{*}** Anti-Static VHC minimizes static charge, increasing the suspension time of particles in the chamber to provide consistent medication delivery for tracheostomy patients.

A Comfortable Fit

The **AeroTrach Plus^{*}** Anti-Static VHC has been designed to quickly connect to typical 15mm adapters. The **Comfort Right Fit^{*}** adapter ensures a firm connection for enhanced aerosol delivery, without the need for excessive torque when connecting or disconnecting the device.

Setting the Gold Standard for Aerosol Delivery

Trudell Medical International is renowned for its extensive research and development of new technology for aerosol delivery. We are proud to continue the tradition of innovation with the **AeroTrach Plus^{*}** Anti-Static VHC, a significant improvement in asthma relief for tracheostomy patients. The **AeroTrach Plus^{*}** Anti-Static VHC is a simple and cost-effective system to deliver medicine where it's needed most – your patient's lungs.



Manufactured by:



Trudell Medical International[®]

725 Third Street, London,
Ontario, Canada, N5V 5G4
Phone: +1-519-455-7060
Fax: +1-519-455-7858

Authorised European Representative:
Trudell Medical International Europe Limited
BioCity Nottingham, Pennyfoot Street
Nottingham, NG1 1GF, UK
Email: customerservice@trudellmed.com
www.trudellmed.com

Canadian Hospital Distribution:

Trudell Medical Marketing Limited



758 Third Street, London,
Ontario, Canada, N5V 5J7
Toll Free: 1-800-265-5494
Email: tmml@tmml.com
www.tmml.com

Ordering Information

Part Number	Description	Case Size
10150591010	AeroTrach Plus[*] Anti-Static VHC	Case of 10

REGISTERED
ISO 13485
QUALITY SYSTEM



ISO CLASS 8
CLEANROOM MANUFACTURING