

the TopChair-S

be free



TopChair

BE FREE



Want to see TopChair-S in action ?

Go to Youtube Channel
www.youtube.com/user/topchair/



TopChair-S

Be **free** !

The power wheelchair is **adapted to all environments : at home, around town or for a countryside stroll**, all without any assistance. It is now possible **to climb obstacles** such as : sidewalks, doorsteps, stairs...

Are you looking for a customizable and powerful wheelchair? **Contact us!**



Road mode :

- **Comfort and security :**
 - The patented technology enables the seated person to remain horizontal while descending or ascending the stairs : **Active Balance System**
 - **Unique shock absorber system :**
 - TFS :** Track Flip Suspension
 - 5 seats available
- **Reliability :**
 - developed and **produced in France** with 2 years warranty
- **Performance :**
 - Speed 9km/h
 - Range 35km (2 batteries 60Ah)
 - Electronic High-tech **R-net** (domotic, Bluetooth, idevices, customization...)

Any specific request?
 We can customize your powerchair!



Smart Climbing mode :

It climbs backward and goes down forwards on straight stairs. While descending or ascending, the wheelchair ensures :

- **Safety :**
 - Seat remains horizontal with automatic : hazard detection (speed, angle obstacle...)
- **Easy to use :**
 - Real time instructions on the screen
- **Certification :**
 - ISO 7176 approval



TopChair has been tested by CERAH
 (Center of study and research for Medical Devices)
and has a CE certification ISO 7176

TopChair-S

Take a closer look

Tilt position

Seats tilt forward 15° to assist user transfer.



Suspension system TFS

Front independent suspension provides a higher level of comfort and stability, indoor and outdoor.



Smart climbing mode

Can go over obstacles with a maximum incline of 33° and 65% and 20 cm maximum height of steps.



Easy to use

- Instructions on the LCD screen (smart climbing mode)
- Swing-away Joystick (road mode)



Colors

- Spicy orange



- Ice white



Technical specifications

LENGTH (without footrest)	115 cm (98cm)
WIDTH	68-70 cm
FRONT WHEELS	230mm /9"
REAR WHEELS	350mm / 14"
SEAT TO FLOOR HEIGHT (without cushion)	45cm (tracks) 56cm (wheels)
WEIGHT OF THE WHEELCHAIR (batteries included)	145kg
MAXIMUM USER WEIGHT	110kg
SEAT WIDTH	Standard 43cm / customizable
SEAT DEPTH	Standard 46cm / customizable
BACKREST WIDTH	Standard 48cm / customizable
BACKREST HEIGHT	Standard 47cm / customizable
ARMSRESTS DISTANCE	35cm – 66cm
ARMREST HEIGHT	2cm – 30cm
BATTERIES TYPE	2x60AH

RANGE	35km
MAXIMUM SPEED	9km/h
TURNING RADIUS (ISO7176-5)	85cm
WIDTH HALF TURN (ISO7176-5)	210cm
MAXIMUM OBSTACLE	Wheels 6cm / tracks 20cm
FRONT SUSPENSIONS	Yes
ELECTRONIC	R-Net
POWER SEAT TILT	30° - 40° backwheels up
POWER BACKREST TILT	Optional
POWER LEGREST	Optional
WARRANTY	2 years
OBSTACLE MAX INCLINE	33° on tracks
MINIMAL TRANSPORT PACKAGE	97cm x 70cm
MINIMAL TRANSPORT HEIGHT	44cm



BE FREE

TopChair SAS

Sumo Park 2 allée Longuetterre
31850 Montrabé / France
Tel : +33 (0)961.227.202
info@topchair.net
www.topchair.fr

FOLLOW US!   

